




Alpha Medical Instruments LLC


BIPOLAR BALLOON CATHETER WITH SHROUDED PIN



Bipolar Pacing Catheter With Shrouded Pin

Series 600

Product Code	Catheter Size	Catheter Body Material	Catheter Body Lumen	Useable Length	Inflated Balloon Diameter	Measurement from Distal Tip Electrode Spacing	Maximum Balloon Inflation
652/1-SH-90P	5 Fr.	PUR	2 L	90 cm	9 mm	1.0 cm	1.0 cc
652/1-SH-110P	5 Fr.	PUR	2 L	110 cm	9 mm	1.0 cm	1.0 cc
652/2-SH-90P	5 Fr.	PUR	2 L	90 cm	9 mm	2.5 cm	1.0 cc
652/2-SH-110P	5 Fr.	PUR	2 L	110 cm	9 mm	2.5 cm	1.0 cc



Description:

The Bipolar Balloon Pacing Catheter is made of extruded radiopaque polymeric tubing with two lumens and two stainless steel electrodes. Electrode spacing is 1.0 or 2.5 cm apart, depending upon which model you choose.

Packaging:

Each catheter is packaged in a sterile peel-pack tray

Catheters Length Marking:

Each catheter has clear, easy to identify increment markings every 10 cm along catheter body