



Alpha Medical Instruments LLC

***DPX-G PEEL AWAY
Introducer Set***



Code: DPX-G

PEEL AWAY INTRODUCER SET



DPX - G Introducer Set

Product Code	Peel Away Dilator Sheath		Syringe	Cannula (Gauge x Length)	Scalpel Size	Guide Wire	
	Size	Length				Length	OD
DPX – 5G	5 Fr	17.3 cm	10 ml	18 Ga. x 7 cm	11	50 cm,	0.038”
DPX – 6G	6 Fr	17.3 cm	10 ml	18 Ga. x 7 cm	11	50 cm	0.038”
DPX – 7G	7 Fr	17.3 cm	10 ml	18 Ga. x 7 cm	11	50 cm	0.038”
DPX – 8G	8 Fr	17.3 cm	10 ml	18 Ga. x 7 cm	11	50 cm	0.038”
DPX – 9G	9 Fr	17.3 cm	10 ml	18 Ga. x 7 cm	11	50 cm	0.038”
DPX – 10G	10 Fr	17.3 cm	10 ml	18 Ga. x 7 cm	11	50 cm	0.038”
DPX – 11G	11 Fr	17.3 cm	10 ml	18 Ga. x 7 cm	11	50 cm	0.038”

Alpha Medical Instruments LLC
 23455 Madero, Suite B
 Mission Viejo, CA 92691
 USA
www.alphamedicalinstruments.com

Tel: 949-460-0944
Fax: 949-460-0904