



**Alpha Medical Instruments LLC**

## ***DPX-T PEEL AWAY Introducer Set***



**Code: DPX-T**

[www.alphamedicalinstruments.com](http://www.alphamedicalinstruments.com)

# ***PEEL AWAY INTRODUCER SET***



# DPX –T PEEL AWAY Introducer Set

Product Code	Peel Away Dilator Sheath		Syringe	Cannula (Gauge x Length)	Scalpel Size	Guide Wire	
	Size	Length				Length	OD
DPX – 5T	5 Fr	19.3 cm	10 ml	18 Ga. x 7 cm	11	50 cm	0.038"
DPX – 6T	6 Fr	19.3 cm	10 ml	18 Ga. x 7 cm	11	50 cm	0.038"
DPX – 7T	7 Fr	19.3 cm	10 ml	18 Ga. x 7 cm	11	50 cm	0.038"
DPX – 8T	8 Fr	19.3 cm	10 ml	18 Ga. x 7 cm	11	50 cm	0.038"
DPX – 9T	9 Fr	19.3 cm	10 ml	18 Ga. x 7 cm	11	50 cm	0.038"
DPX – 10T	10 Fr	19.3 cm	10 ml	18 Ga. x 7 cm	11	50 cm	0.038"
DPX – 11T	11 Fr	19.3 cm	10 ml	18 Ga. x 7 cm	11	50 cm	0.038"

**Alpha Medical Instruments LLC**  
 23455 Madero, Suite B  
 Mission Viejo, CA 92691  
 USA  
[www.alphamedicalinstruments.com](http://www.alphamedicalinstruments.com)

**Tel: 949-460-0944**  
**Fax: 949-460-0904**